
Carboxymethylcellulose sodium salt, Medium viscosity (CMC)

RM10844

Product Identifier

| | | |
|------------|---|------------|
| CAS No. | : | 9004-32-4 |
| HS Code | : | 3912 31 00 |
| Storage | : | Below 30°C |
| Shelf life | : | 4 years |

Technical Specification

| | | |
|---------------------------------|---|-------------------------------------|
| Appearance | : | White to yellow to light tan powder |
| Solubility | : | 10 mg soluble in 1 mL of water |
| pH (1% in water at 25°C) | : | 6.50 - 8.00 |
| FTIR | : | Matches with the standard pattern |
| Loss on drying (at 110°C, 2 hr) | : | <= 15.00% |
| Viscosity (2% in water at 25°C) | : | 400 - 800 cps |

Safety Information

| | | |
|---------------|---|---------------------|
| UN No. | : | Not dangerous goods |
| Class | : | - |
| Packing Group | : | - |
| RTECS | : | FJ5950000 |
| WGK | : | 1 |