

Product Specification

Product Name:
Benzyl chloroformate - (stabilised) for synthesis

Product Number: 8.14688
CAS Number: 501-53-1
Formula: C₈H₇ClO₂
Formula Weight: 170,59 g/mol

TEST	Specification
Assay (Acidimetric)	≥ 96.0 %
Density (D ₂₀ /4)	1.205 - 1.215
Infrared Spectrum	Passes Test

Remarks:

The material contains 0.1% anhydrous sodium carbonate as stabiliser.

Specification: PRE.0.ZQ5.10000025783

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

