



Specification

1.06237.0000 N-(1-Naphthyl)ethylenediamine dihydrochloride GR for analysis ACS

	Specification	
Assay (argentometric, calc. on anhydrous substance)	≥ 97.0	%
Identity (IR-spectrum)	conforms	
Appearance of solution (20 g/l, water)	clear, dissolution is complete	
1-Naphthylamine (HPLC)	≤ 0.1	%
2-Naphthylamine (HPLC)	≤ 0.01	%
Water (according to Karl Fischer)	≤ 5.0	%
Sensitivity to sulfanilamide	conforms	

Dr. Sebastian Lips

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.