



Specification

1.08114.0000 Tetrahydrofuran EMPLURA®

	Specification	
Purity (GC)	≥ 99.0	%
Identity (IR)	conforms	
Free acid (as CH ₃ COOH)	≤ 0.002	%
Density (d 20 °C/ 4 °C)	0.887 - 0.889	
Peroxide (as H ₂ O ₂)	≤ 0.01	%
Water	≤ 0.1	%

Stabilized with 2,6-Di-tert-butyl-4-methylphenol (BHT)

Jeannette David

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.