
Pararosaniline hydrochloride, Hi-AR™/ACS

RM3059

Product Identifier

CAS No.	:	569-61-9
EC No.	:	209-321-2
C.I. No.	:	42500
Molecular Formula	:	C ₁₉ H ₁₇ N ₃ .HCl
Molecular Weight	:	323.82
HS Code	:	3204 13 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Metallic green powder.
Solubility	:	Soluble in methanol.
Assay (UV)	:	min. 99.00%

Safety Information

Hazard Pictogram(s) :



Signal Word	:	Danger
Hazard Statement(s)	:	H350
Precautionary Statement(s)	:	P201- P308 + P313
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	CX9850100
WGK	:	3