

RM7625 Xylitol

Product Number Packing

RM7625 : 100G
RM7625 : 25G
RM7625 : 5G

Product Information

Product Code : RM7625
Product Name : Xylitol
Synonym : Xylite
Molecular Formula : C₅H₁₂O₅
Molecular Weight : 152.15
CAS No. : 87-99-0
Shelf Life : 4 years

Technical Specification

Appearance : Colourless to white, hygroscopic crystals or crystalline powder.
Solubility : Soluble in water.
Melting range : 92 - 96°C
Water (K.F.) : <= 0.5%
Assay (HPLC) : min. 98.00 %

Risk and Safety Information

RTECS : ZF0800000
Storage Temperature(°C) : Store below 30°C

Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods