



Lactophenol Cotton Blue

S016

Intended use

Lactophenol Cotton Blue is used for staining solution for fungi.

Composition**

Ingredients

Phenol crystals	20.0 gm
Cotton blue	0.05 gm
Lactic acid	20.0 ml
Glycerol	20.0 ml
Distilled water	20.0 ml

**Formula adjusted, standardized to suit performance parameters

Directions

1. Place a drop of Lactophenol Cotton Blue reagent on a clean and dry slide. The stain imparts a blue colouration on hyphae.
2. By using a nichrome inoculating wire, carefully tease the fungal culture into a thin preparation.
3. Place a coverslip on the preparation. Wait for about 5 minutes.
4. Observe first under microscope with low power for screening in low intensity.

Principle And Interpretation

Fungus are eukaryotic organism and they are classified into two main groups that is yeast and molds. Its cell wall is made up of chitin. Fungal cells have both macroscopic as well as microscopic structure. Lactophenol along with Cotton Blue reagent is used for staining as well as for wet mounting of fungi. Wet mount preparation is the most widely used method of staining and observing fungi and is simple to prepare. The preparation has three components: Lactic acid preserves the fungal structure and clears the tissue while phenol acts as a disinfectant and cotton blue imparts blue colouration to the fungal spores and hyphae.

Type of specimen

Primarily with pure cultures ; Clinical samples - skin, hair and nail tissue (Certain specimens may be examined directly using this stain, others may require processing).

Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines. After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic Use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and

culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. Lactophenol Cotton Blue is useful in the recognition and presumptive identification of fungi. Additional characteristics including colony morphology and biochemical tests should be used where appropriate for final identification
2. Disruption of the fragile fungal architecture during sampling may occur.
3. Wet mount staining do not allow observation of early stage differentiation of the fungus.
4. Wet mount can not be stored over longer period of time.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature

Quality Control

- **Appearance** : Ink blue coloured solution.
- **Clarity** : Clear, without any insoluble particles.
- **Microscopic Examination** : Fungal staining is carried out. Fungal Spores and hyphae are observed under microscope using high power (40X) objective lens after staining with Lactophenol cotton blue.
- **Results** : Fungal spores and hyphae : pale to dark blue.

Storage and Shelf Life

Store between 10-30°C in tightly closed container and away from bright light. Use before expiry date on label. On opening, product should be properly stored in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques.

Reference

1. Generalmicroscience.com/microbiollaboratorytechniques/ staining fungus using lactophenol cotton blue.
2. Koneman, E.W., et al. Color Atlas and Textbook of Diagnostic Microbiology, J.B. Lippincott Company, Philadelphia, PA.
3. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
4. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.



Storage temperature



Do not use if package is damaged



In vitro diagnostic medical device



CE Marking



HiMedia Laboratories Pvt Limited
C-40,21/Y, MIDC, Wagle Ind Area,
Thane(W)-400604,Maharashtra,India



CEpartner4U,ESDOORNLAAN 13,3951
DB MAARN,The Netherlands,
www.cepartner4u.eu

Revision 01/2022

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.