

Product Specification

Product Name:

Potassium chromate - for analysis EMSURE® ACS,Reag. Ph Eur

Product Number:

1.04952

CAS Number:

7789-00-6

Formula:

K₂CrO₄

Formula Weight:

194,2 g/mol

TEST	Specification
Assay (Iodometric)	≥ 99.5 %
Insoluble Matter (In Water)	≤ 0.005 %
pH (5 %; Water)	9.0 - 9.8
Calcium (Ca)	≤ 0.005 %
Sodium (Na)	≤ 0.02 %
Lead (Pb)	≤ 0.005 %
Chloride (Cl)	≤ 0.001 %
Sulfate (SO ₄)	≤ 0.01 %

Specification: PRE.0.ZQ5.10000026983

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

