



Specification

1.06268.0000 Sodium acetate anhydrous for analysis EMSURE® ACS,Reag. Ph Eur

Spec. Values

Assay (perchloric acid titration)	≥ 99.0	%
Identity	passes test	
Appearance of solution	passes test	
Insoluble matter	≤ 0.01	%
pH-value (5 %; water, 25 °C)	7.0 - 9.2	
Chloride (Cl)	≤ 0.002	%
Phosphate (PO ₄)	≤ 0.001	%
Sulfate (SO ₄)	≤ 0.003	%
Al (Aluminium)	≤ 0.001	%
Ca (Calcium)	≤ 0.005	%
Cu (Copper)	≤ 0.0003	%
Fe (Iron)	≤ 0.001	%
K (Potassium)	≤ 0.05	%
Mg (Magnesium)	≤ 0.002	%
Heavy metals (ACS)	≤ 0.001	%
Loss on drying (120 °C)	≤ 1.0	%

Corresponds to ACS,Reag. Ph Eur

Kristin Starke

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.