

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
(3-Aminopropyl)triethoxysilane - $\geq 98.0\%$

Product Number: **741442**

CAS Number: 919-30-2



Formula: C₉H₂₃NO₃Si
Formula Weight: 221.37 g/mol

TEST

Specification

Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared Spectrum	Conforms to Structure
Titration with HCl	98.0 - 102.0 %
Color Test	≤ 25 APHA
Purity (GC)	≥ 98.0 %

Specification: PRD.3.ZQ5.10000020462

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.