

**GRM574 Naphthalene, Hi-LR™****Product Number Packing**

GRM574 : 500G

**Product Information**

Product Code : GRM574  
Product Name : Naphthalene, Hi-LR™  
Molecular Formula : C<sub>10</sub>H<sub>8</sub>  
Molecular Weight : 128.17  
CAS No. : 91-20-3  
EC No. : 202-049-5  
HS Code : 2902 90 00  
EC Index No. : 601-052-00-2  
Shelf Life : 4 years

**Technical Specification**

Appearance : White to almost white crystals or powder or flakes or balls  
Solubility : 33.3 mg soluble in 1 mL of ethanol  
FTIR (KBr disc) : Matches with the standard pattern  
Melting range : 79 - 82°C  
Assay (GC) : min. 99.00 %

**GHS Safety Information**

Hazard Statement(s) : H302-H351-H410  
Precautionary Statement(s) : P273-P281-P501  
Signal Word : Warning  
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 22-40-50/53  
S-Phrase(s) : 36/37-46-60-61  
WGK : 3  
RTECS : QJ0525000  
Flash Point(°F) : 176 °F  
Flash Point(°C) : 80°C  
Storage Temperature(°C) : Store below 30°C  
Hazard Symbol(s)

**Transport Information**

UN No. : 1334



## Product Information

Class	:	4.1
Packaging Group	:	3
Marine Pollutant	:	Yes
ADR/RID	:	1334 4.1/PG 3
IMDG	:	1334 4.1/PG 3
IATA	:	1334 4.1/PG 3

A-406, Bhaveshwar Plaza, LBS Marg, Mumbai-400 086, India. Phone:+91-022-4095 1919 Fax:+91-022-4095 1920 Gram: STERILITY.  
email:info@himedialabs.com. Visit us at our website: [www.himedialabs.com](http://www.himedialabs.com)

---

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein