

03.1031.00 - 03.1035.00



**ANELLI PER SOSTEGNO IN FERRO VERNICIATO, CON MORSETTO A VITE**

Apertura morsetto 18 mm, interasse tra morsetto e anello 120 mm. Sono consigliati come supporti per imbuto separatori sferici.

**PAINTED IRON RINGS FOR STAND WITH SCREW CLAMP**

Clamp opening 18 mm, distance between clamp and ring 120 mm. We suggest to use them as supports for spherical separating funnels.

Cod.	Diametro   Diameter	Lunghezza   Length
03.1031.00	50 mm	
03.1031.01	70 mm	120 mm
03.1032.00	85 mm	
03.1033.00	115 mm	125 mm
03.1034.00	135 mm	
03.1035.00	150 mm	140 mm

03.1036.00 - 03.1040.00



**ANELLI PER SOSTEGNO IN ACCIAIO INOX, CON MORSETTO A VITE**

Apertura morsetto 18 mm, interasse tra morsetto e anello 120 mm. Sono consigliati come supporti per imbuto separatori sferici.

**STAINLESS STEEL RINGS FOR STAND, WITH SCREW CLAMP**

Clamp opening 18 mm, distance between clamp and ring 120 mm. We suggest to use them as supports for spherical separating funnels.

Cod.	Diametro   Diameter	Lunghezza   Length
03.1036.00	50 mm	
03.1036.01	70 mm	120 mm
03.1037.00	85 mm	
03.1038.00	115 mm	125 mm
03.1039.00	135 mm	
03.1040.00	150 mm	140 mm

03.1042.00



**ANELLI OPPOSTI PER SOSTEGNO**

In ottone cromato, con morsetto centrale per aste fino a 10 mm di diametro, adatti per filtrare.

**OPPOSITE RING FOR SUPPORT**

Made of chromate brass, with central screw clamp for rod up to 10mm diameter, suitable for filtering.

Diametro interno | Internal diameter 45 mm

03.1044.00



**ANELLI DOPPI PER SOSTEGNO**

In ottone cromato, con due anelli posizionati a 120°, con morsetto centrale per aste fino a 10 mm di diametro, adatti per filtrare.

**DOUBLE RINGS FOR STAND**

Made of chromate brass, with two rings positioneted at 120°, with central screw clamp for rod up to 10mm diameter, suitable for filtering.

Diametro interno | Internal diameter 45 mm

03.1169.00



**MORSETTO PER AGITATORI AD ASTA**

In lega anticorrosiva, apertura 6 - 18 mm.

**DOUBLE CLAMP FOR STIRRER ROD**

Made of antiacid and anticorrosive light alloy, opening 8 - 18 mm.

03.1175.00



**MORSETTO A CROCE UNIVERSALE**

In alluminio verniciato, apertura 16 mm.

**UNIVERSAL CROSS CLAMP**

Made of painted aluminum, opening 16 mm.