



Specification

1.07838.0000 Antimony(III) chloride for analysis EMSURE® ACS

	Specification	
Assay (bromometric)	≥ 99.0	%
Matter insoluble in chloroform	≤ 0.02	%
Substances insoluble in hydrochloric acid	≤ 0.005	%
Sulfate (SO ₄)	≤ 0.005	%
As (Arsenic)	≤ 0.005	%
Ca (Calcium)	≤ 0.005	%
Cu (Copper)	≤ 0.0005	%
Fe (Iron)	≤ 0.002	%
K (Potassium)	≤ 0.01	%
Na (Sodium)	≤ 0.005	%
Pb (Lead)	≤ 0.003	%

Dr. Dimitrij Ryvlin

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.