

Oryzalin, Plant Culture Tested

PCT1559

Product Identifier

CAS No.	:	19044-88-3
EC No.	:	242-777-0
Molecular Formula	:	C ₁₂ H ₁₈ N ₄ O ₆ S
Molecular Weight	:	346.36
HS Code	:	2935 90 90
Storage	:	Below 30°C
Shelf life	:	2 years

Technical Specification

Appearance	:	Yellow to orange crystals or powder or solid.
Solubility	:	1 mg soluble in 1 mL of methanol
Assay (HPLC)	:	min.95.00%
Plant Tissue Culture Test	:	Passes

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H317- H351- H373- H410
Precautionary Statement(s)	:	P201- P260- P273- P280- P308 + P313- P314- P501
UN No.	:	3077
Class	:	9
Packing Group	:	III
RTECS	:	WO9350000
WGK	:	1