

## FEATURES

- Shaping shell is made of cold rolled steel sheet with surface spraying treatment
- Nichrome heating wire, faster warming up time, long service life
- Easy operation, and wide application
- Precise PID temperature control, meets high accuracy experiment needs

## PRODUCT INFORMATION

- Product Name: Several Row Heating Mantle with Soxhlet
- Model Number: SASE-250-4, SASE-500-4, SASE-1000-4  
SASE-250-6, SASE-500-6, SASE-1000-6  
SDSE-250-3, SDSE-500-3, SDSE-500-3  
SDSE-250-6, SDSE-500-6, SDSE-500-6
- Brand: B-ONE
- Product Category: Laboratory instrument – Heating mantle

## APPLICATIONS

- Chemical Reactions
- Distillation and Evaporation
- Extraction
- Synthesis
- Heating

## SPECIFICATIONS

### SEVERAL ROW ANALOG HEATING MANTLE

MODEL	CAPACITY	MAX TEMPERATURE	VOLTAGE	POWER
SASE-250-4	4 x 250ml	450 °C	220VAC/50Hz	4 x 150W
SASE-500-4	4 x 500ml			4 x 250W
SASE-1000-4	4 x 1000ml			4 x 350W
SASE-250-6	6 x 250ml			6 x 150W
SASE-500-6	6 x 500ml			6 x 250W
SASE-1000-6	6 x 1000ml			6 x 350W

### SEVERAL ROW DIGITAL HEATING MANTLE

MODEL	CAPACITY	MAX TEMPERATURE	VOLTAGE	POWER
SDSE-250-3	3 x 250ml	450 °C	220VAC/50Hz	4 x 150W
SDSE-500-3	3 x 500ml			4 x 250W
SDSE-1000-3	3 x 1000ml			4 x 350W
SDSE-250-6	6 x 250ml			6 x 150W
SDSE-500-6	6 x 500ml			6 x 250W
SDSE-1000-6	6 x 1000ml			6 x 350W

## INCLUDES

- Soxhlet
- Extraction glassware
- Clamp
- Statif pole
- Bosehead

## WARRANTY

- One year (parts & service)
- Warranty does not apply to damages caused by user error, misuse, or failure to follow operating instructions.



Several Analog Heating Mantle with Soxhlet



Several Digital Heating Mantle with Soxhlet