


Product Specification

Product Name:
N,O-Bis(trimethylsilyl)trifluoroacetamide - ≥99%

Product Number: **155195**
CAS Number: 25561-30-2



Formula: C₈H₁₈F₃NOSi₂
Formula Weight: 257,40 g/mol

TEST	Specification
Appearance (Color)	Colorless to Almost Colorless
Appearance (Form)	Liquid
Color Test	≤ 50 APHA
Purity (GC)	≥ 99,0 %
Infrared Spectrum	Conforms to Structure

Specification: PRE.0.ZQ5.10000021368

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

